

KOHLA

EN


CLIMBING SKIN

SKINFINDER

HOW TO CHOSE THE RIGHT SKIN

ALPINIST	VERTICAL	VERTICAL WMN	FREERIDE FREERIDE AIR	PEAK
100% Mohair	Mixmohair	Mixmohair	Mixmohair	Mixmohair
Fiberseal	Fiberseal	Fiberseal	Fiberseal	Baumwolle
Smart Glue	Smart Glue / Vacuum	Smart Glue	Smart Glue	Hotmelt
130 mm width	120 / 130 mm width	135 mm width	135 mm width	120 /135 mm width

1. Measure actual ski length
2. Measure the width of the front ski tip
3. Select the right skin in the table
(For the right skin length: ski length -1cm)



ski length (cm)	skin length (cm)	width		
		up to 120 mm	up to 130 mm	up to 135 mm
142 - 148	142	Vertical / Vertical wmn / Peak		
149 - 155	149	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
156 - 162	156	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
163 - 169	163	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
170 - 176	170	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
177 - 183	177	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
184 - 190	184	Vertical / Vertical wmn / Peak	Alpinist	Freeride / Freeride AIR / Peak
191 - 197	191		Alpinist	Freeride / Freeride AIR / Peak

FREERIDE AIR KIDS	ski length (cm)	skin length (cm)
Mixmohair	125 - 134 cm	125
Fiberseal	135 - 144 cm	135
Smart Glue	145 - 154 cm	145
120 mm width	155 - 164 cm	155

SCAN ME



KOHLA SERVICE
VIDEOS